# SLStudio User Manual

Version 1.0, Revised 24/04/2021


By: Lam Teng Foong

### Required Software

- QT Creator 5.9.5
- Point Cloud Library (PCL) 1.8
- OpenCV OpenCV 3.2.0
- ROS with the following packages
  - Versavis (devel/sl\_sensor branch from <u>forked repo</u>)
  - flir\_camera\_driver (devel/sl\_sensor branch from forked repo)

# Running SLStudio

- Start Versavis to begin image acquisition
  - If using hardware trigger:
    - In versavis\_configuration.h, change the camera trigger frequency based on the type of pattern used (PSP+TPU or 2P1 with TPU). Header file has comments to indicate with lines of code to comment/uncomment
  - Upload versavis.ino onto the board
  - In a terminal, run the code
    - roslaunch versavis sls\_run\_versavis.launch
- Launch SLStudio
  - In a separate terminal (need to be launch for a terminal to run properly with ROS), run
    - qtcreator
  - Open the SLStudio project in qtreator
  - Press Run (Large green triangle)



#### Preferences


Pattern to be projected. Currently only supports:

- 2x3 Pattern Phase Shift
- 2+1 Pattern with TPU

Project vertical and/or horizontal patterns. Note 2+1 TPU only supports vertical or horizontal, but not both together

Check this because Lightcrafter has a diamond pixel arrangement

Shutter option is currently not in use, exposure time is managed by flir\_camera\_driver ROS package and is set to 16.667ms



#### Calibration

SL Calibration All captured ✓ Sequence 8 ✓ Sequence 9 sequences, untick low ✓ Sequence 10 ✓ Sequence 11 quality sequences to √ Sequence 12 √ Sequence 13 remove them from the √ Sequence 14 Screen √ Sequence 15 calibration process √ Sequence 16 showing ✓ Sequence 17 Live View / √ Sequence 18 Sequence 19 **Processed** ✓ Sequence 20 ✓ Sequence 21 Calibration ✓ Sequence 22 Calibration board √ Sequence 23 images grid size Row and columns Size (mm): 8 of calibration 10 Rows: board (# of grid 10 Cols: intersections) Calibrate Press Calibrate to Cancel Save estimate sensor Snap Re-Calibrate intrinsics/extrinsics Press Save once you are satisfied with the Snap to take a sequence results. Results will be saved in the file of calibration images calibration.xml

### Post-Calibration Steps

- 1. Once calibration has been done, press 'Generate Patterns' in the main menu
- 2. Upload the generated patterns onto the Lightcrafter projector using the steps in <u>pattern upload instructions.pdf</u>

# Linearity Testing

- 1. Upload the DLPR350PROM\_v4.1.0\_with\_cali\_linearity\_test.bin firmware onto the projector using the <u>lightcrafter 4500 gui</u>
- 2. Press 'Linearity Test' in the main menu
- 3. Projector will project fully illuminated images of increasing intensity and the camera will take pictures of them
- 4. Use the 'utiliy/linearity\_test.m' script in <a href="mailto:phase shifting profilometry matlab">phase shifting profilometry matlab</a> to evaluate the captured images